



Introducing

Quantum™ 2

COBLATION®

For Intelligent Arthroscopy


The Quantum 2 Difference

ArthroCare® invented Coblation® for arthroscopic surgery, and no one understands the science or technology better. Our engineers have taken everything we've learned over the past 15 years and more than 5 million procedures to create the Quantum 2 System, our revolutionary new Coblation platform. Quantum Coblation takes arthroscopic surgery to an entirely new level by giving you the control you need to handle whatever you find, in every joint you treat.

How Quantum Coblation Works

Low Temperature Plasma Excision

Quantum Coblation uses a controlled, non-heat driven process in which bipolar radiofrequency (bRF) energy excites the electrolytes in a conductive medium, usually saline solution, to create a precisely focused and charged plasma gas. The energized particles in the plasma have sufficient energy to break the molecular bonds within tissue, causing tissue to dissolve at relatively low temperatures (between 40°C and 70°C). Because radiofrequency current does not pass directly through tissue during Coblation, tissue heating is minimal. The result is volumetric removal of the target tissue with minimal damage to surrounding healthy tissue.



During Coblation, high voltage is applied to the conductive medium (usually saline) located between the electrode and the tissue. The high voltage converts the fluid into an ionized vapor layer (plasma).



The plasma contains excited particles that accelerate toward the tissue and break the tissue's molecular bonds. At relatively low temperatures (between 40°C and 70°C), the tissue disintegrates molecularly (not explosively) into simple molecules, such as small hydrocarbons and gases (oxygen, nitrogen, carbon dioxide, etc.) that are removed from the site with suction, resulting in effective tissue removal.

The Quantum 2 System is the most dynamic controller we've created. When paired with our Multi-Electrode ArthroWand™ devices, the controller produces a large and responsive plasma layer for fast ablation speeds. The Quantum 2 System offers a unique set of features designed to make your procedures go smoothly.

Efficient Ablation, Intelligent Features



Multi-Electrode Technology.™ The ArthroWand portfolio offers a full range of Multi-Electrode wands. These wands form a stable plasma layer, which leads to fast ablation and responsive performance.



Hand Control Capability. The Quantum 2 System supports a line of wands with Integrated Finger Switches (IFS).




Temperature Feedback. The Ambient Wand Collection provides accurate, real-time temperature measurements of the operative environment by using visual and audio indicators on the Quantum 2 Controller.



Wireless Footswitch. The Quantum 2 System is compatible with a wireless footswitch accessory which eliminates the tangle of controller wires for unsurpassed convenience and ease of use.



Adjustable Coagulation. The new enhanced coagulation mode can improve visualization when tough bleeders are encountered.



Automatic Scope Saver. The Quantum 2 System automatically pauses radiofrequency output when wands are in close proximity to or in contact with metal, and resumes radiofrequency output when the wands have returned to a proper distance.




Integrated TOPAZ™ Timer. When using TOPAZ Wands, the Quantum 2 Controller defaults to optimal power and timer settings, automatically delivering the precise dose.

Comprehensive Wand Technologies

The flexible and versatile Quantum 2 System was designed to meet all of your sports medicine needs. Our extensive line of ArthroWand devices offers the widest variety of angles and electrode configurations available. Whether you're debriding articular cartilage in the knee, performing hip arthroscopy, or treating a tendon, Quantum Coblation and our innovative ArthroWand devices are the best choice for the most demanding procedures.

The Right Tools for Every Indication



Hip Arthroscopy

SideWinder™ Blade, the first articulating bipolar hip Wand, incorporates a rotating shaft that is deflectable up to 100° enabling unmatched access in the joint space. MultiVac™ 50 XL and TurboVac™ 90 XL offer the length desirable for hip ablation. The hip Wand collection provides rapid removal of synovial tissue and necessary sculpting precision for acetabular labral tissue.



Chondroplasty

The Paragon T2™ Wand is based on 6 years of published, scientific research on the application of bipolar radiofrequency energy on articular cartilage. Temperature Technology, or T2, provides a real-time visual indicator when the temperature exceeds the recommended 50° C.



Tendon and Fascia Treatment

TOPAZ Microdebrider™ devices are designed for the minimally invasive treatment of tendons and fascia, helping individuals return to sports, work, and life. The small tip helps to preserve tissue structure.

The ArthroWand Collection

For 15 years, ArthroCare Sports Medicine has advanced arthroscopic Wand design. Our portfolio of ArthroWand devices includes tools for shoulder, knee, hip, ankle, elbow, wrist, and tendon procedures.



StarVac™ 90



Super TurboVac™ 90



UltraVac™ 90



TurboVac™ 90



Paragon T2



Super MultiVac™ 50



CoVac™ 50



TriStar™ 50




CoVator™ 20



SideWinder



Bevel 45




CAPSure™ 30



Microblator™ 30



TOPAZ MicroDebrider



LoPro™ 90



Eliminator™ 90



Dome 60














Saber™ 30

ArthroWand Selection Chart



RECOMMENDED SETTINGS

	Catalog No.		Wand	Shaft Size	Angle	5.5mm ² Cannulas	KNEE							
	Non-ICW ¹	ICW ¹					VAR	6-8	7-9	6-7	5-7	6-9	4	1-2
Suction 		ASC4251-01	StarVac™	3.75mm	90°	X			■			■		
		ASHA4250-01	Ambient™ Super TurboVac™ 90 IFS ⁴	3.75mm	90°				■			■		
		ASH4250-01	Super TurboVac 90 IFS ³	3.75mm	90°	X			■			■		
		ASC4250-01	Super TurboVac 90	3.75mm	90°	X			■			■		
		ASHA4830-01	Ambient Super MultiVac™ 50 IFS ⁴	3.75mm	50°		7-9		■			■		
		ASH4830-01	Super MultiVac 50 IFS ³	3.75mm	50°		7-9		■			■		
		ASC4830-01	Super MultiVac 50	3.75mm	50°		7-9		■			■		
		AS1335-01	ASC1335-01	TurboVac™ 90	3.75mm	90°	X			■			■	
		AS1337-01		TurboVac 90 HP	3.75mm	90°								
		AS2530-01	ASC2530-01	CoVac™ 50	3.0mm	50°		6-8					■	
	A S4630-01	ASC4630-01	TriStar™ 50	3.0mm	50°	X	6-8			■		■		
		ASC5000-01	UltraVac™ 90	4.25mm	90°				■			■		
Hip 		ASC1336-01	TurboVac 90 XL	3.75mm	90°	X								
		ASC4730-01	MultiVac 50 XL	3.75mm	50°	X								
		AC2340-01	SideWinder™ Blade	2.7mm	55°	X								
Right Angle 	A1325-01		Right Angle 90	2.5mm	90°						■			
	A1335-01	AC1340-01	Right Angle 90	3.5mm	90°				■			■		
	A1336-01	AC1336-01	LoPro™ 90	3.6mm	90°	X			■			■		
	A1345-01	AC1345-01	Eliminator™ 90	4.5mm	90°	X			■			■		
Tendon 		ACH4040-01	TOPAZ™ MicroDebrider IFS ³	0.8mm	0°								■	
		AC4040-01	TOPAZ MicroDebrider	0.8mm	0°								■	
		ACH4045-01	TOPAZ XL MicroDebrider IFS ³	1.27mm	15°	X								
		AC4045-01	TOPAZ XL MicroDebrider	1.27mm	15°	X								
Foot & Ankle 		Q6000-01	TOPAZ MicroDebrider	0.8mm	0°								■	
		Q6002-01	TOPAZ EPF MicroDebrider	0.8mm	45°									
		Q6003-01	TOPAZ EPF Cannula System											
Cutting 	A4330-01	AC4330-01	Saber™ 30	3.0mm	30°				■					
		AC4340-01	CoVator™ 20	3.3mm	20°	X			■					
Bevel 	A2430-01	AC2430-01	Bevel 45	3.0mm	45°	X	6-8		■		■			
	A2530-01		Bevel 60	3.0mm	60°		6-8		■			■		
	A2630-01		Bevel 30	3.0mm	30°	X	6-8							
Dome 	A3525-01	AC3525-01	Dome 60	2.5mm	60°		6-8							
	A3625-01		Dome 30	2.5mm	30°	X	6-8							
Meniscus & Articular Cartilage 		AC5520-01	Razor2.5™ 45	2.5mm	45°	X	4-9							
		AC5531-01	Paragon T2™	2.3mm	15°	X					■			
Coagulation 	A1720-01		Micro CAPS™ ↗	2.0mm	0°	X								
		AC1830-01	CAPSure 30 ↗ ⊙	3.0mm	30°	X							■	
Small Joint 		AC2823-01	Short Bevel 35	2.3mm	35°									
		AC4050-01	Microblator™ 30	1.4mm	30°									

Meniscectomy	Lateral Release	ACL/PCL Debridement (Reticuloplasty)	Articular Cartilage Debridement (Chondroplasty)	Synovectomy/Plica Removal	Tendon Debridement	Ligament/Tendon Coagulation
--------------	-----------------	--------------------------------------	---	---------------------------	--------------------	-----------------------------

Count on ArthroCare For Leadership and Innovation

For 15 years, ArthroCare has been the leading innovator in radiofrequency solutions for arthroscopic surgery. No one understands the technology better. No one provides a broader range of solutions. No one has an organization so dedicated to these systems. Let us show you the difference Quantum Coblation can make for your practice.

Visit us online at
quantumcoblation.com



International

ArthroCare Europe AB
Skeppsbron 2
111 30 Stockholm
Sweden
Tel: +46 8 546 172 00
Fax: +46 8 546 172 39
sportsmedicine@arthrocare.se

Asia Pacific

ArthroCare (Hong Kong) Limited
Unit 3328, 33rd Floor,
China Merchants Tower
Shun Tak Centre
168 Connaught Road Central
Hong Kong
Tel: +852 3793 3430/32
Fax: +852 3793 3368

Australia

ArthroCare (Australasia) Pty Ltd
PO Box 2186
Normanhurst NSW 2076
Australia
Tel: +61 2 9479 9200
Fax: +61 2 9479 9299
service@arthrocare.net.au

Austria

ArthroCare Austria AG
Kaufmannstrasse 16
6020 Innsbruck
Austria
Tel: +43 512 36 15 38-00
Fax: +43 512 36 15 38-60
info@arthrocare.at

France

ArthroCare France SARL
9, rue Noiroit
70000 Vesoul
France
Tel: +33 3 84 76 81 38
Fax: +33 3 84 76 74 63
info@arthrocare.fr

Germany

ArthroCare (Deutschland) AG
Berghauserstrasse 62
42859 Remscheid
Germany
Tel: +49 2191 93342 0
Fax: +49 2191 93342 49
info@arthrocare.de

Scandinavia

ArthroCare Danmark ApS
Science Center Skejby
Brendstrupgårdsvej 21
DK-8200 Århus N.
Denmark
Tel: +45 7023 2850
Fax: +45 7023 2851
info@arthrocare.dk

Singapore

ArthroCare Singapore Pte Ltd
120 Lower Delta Road
Cendex Centre #12-11
Singapore 169208
Tel: +65 6276 9022
Fax: +65 6278 7622
customer.support@arthrocare.com

Switzerland

ArthroCare Europe AB,
Zuchwil Branch (Switzerland)
Gewerbstrasse 2
CH-4528 Zuchwil
Switzerland
Tel: +41 32 686 88 99
Fax: +41 32 686 88 90
info@arthrocare.ch

UK/Ireland

ArthroCare UK Ltd
St James Business Park
9, Grimbold Crag Court
Knaresborough
North Yorkshire HG5 8QB
United Kingdom
Tel: +44 1423 888806
Fax: +44 1423 888816
info@arthrocareuk.com

USA

ArthroCare Sports Medicine
680 Vaqueros Avenue
Sunnyvale, CA 94085-3523
United States
Tel: +1 408-736-0224

Not for distribution in the United States.

© 2010 ArthroCare Corporation. All rights reserved. All trademarks, service marks, and logos listed herein are registered and/or unregistered trademarks of ArthroCare Corporation, its affiliates, subsidiaries, or a third party who has licensed its trademarks to ArthroCare. ArthroCare ENT, ArthroCare Sports Medicine, and ArthroCare Spine are affiliates of ArthroCare Medical Corporation.