



LubriPen

*KLS Martin presents the new LubriPen
for precise lubrication of surgical instruments
and medical devices*



- instrument oil, suitable for steam and hot air sterilization
- precise dosed without spraying
- physiologically safe
- odourless
- colourless
- steam permeable
- temperature stability up to 260°C
- silicone-free
- leaves no residuals on the instrument surfaces
- provides an excellent corrosion and oxidation protection
- meets the requirements of DAB, LFGB, FDA and USDA (H1)
- 7 mm stainless steel tip, inner diameter 0.8 mm
- pen contents 12 ml or 1,000 drops

**Available on stock.
Just order 55-997-01-04**



03.16 · 90-102-02-04