

# Aesculap® AdTec® Noir®

Endoscopic Noir® Scissors



Aesculap Endoscopic Technology

# Aesculap Endoscopic Technology

## AdTec® Noir®



Noir® stands for  
„No Irritating Reflections“

Noir® steht für  
„No Irritating Reflections“

### Glare-free operation

Black coating of the scissor blades against light reflections – for irritation-free surgery.

### No darkening of camera images

The black scissor blades prevent unintentional darkening of the video image by the camera (Auto Shutter function).

### Reduced tissue sticking

Smoother Noir® surface coating for reduced coagulation tissue sticking.

### Long service life

with the well-tried Aesculap carbide inserts.

### Blendfreies Arbeiten

durch die schwarz beschichteten Scherenblätter werden Lichtreflexionen verhindert. Keine Irritationen des Operateurs.

### Kein Abdunkeln des Kamerabildes

die schwarzen Scherenblätter verhindern, dass die Kamera das Videobild ungewollt abdunkelt (Auto Shutterfunktion).

### Reduziertes Anhaften von Gewebe

durch glattere Noir®-Oberflächenbeschichtung. Gewebeanhaftungen bei der Koagulation werden vermindert.

### Lange Nutzungsdauer

durch die bewährten Aesculap Hartmetalleinlagen.



**PO090R**

Noir® METZENBAUM

scissors, serrated, tips curved to the left, with carbide inserts  
 Schere, gezahnt, Blätter nach links gebogen, mit Hartmetalleinlage



1/2

Single Parts / Einzelteile

P0790R

PM973R

P0958R



**PO091R**

Noir® Mini-METZENBAUM

scissors, serrated, tips curved to the left, with carbide inserts  
 Schere, gezahnt, Blätter nach links gebogen, mit Hartmetalleinlage



1/2

Single Parts / Einzelteile

P0791R

PM973R

P0958R



Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 (0) 74 61 95-0 | Fax +49 (0) 74 61 95-26 00 | [www.aesculap.com](http://www.aesculap.com)

Aesculap – a B. Braun company

The main product trademark 'Aesculap' and the product trademarks 'AdTec' and 'Noir' are registered trademarks of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.

Brochure No. C23711

0911/1.0/3