



D > SHAVERBLADES

GB > SHAVERBLADES

F > SHAVERBLADES

E > SHAVERBLADES

AGGRESSIVE FULL RADIUS RESECTOR
innen gezahnt - außen ungezahnt
inner blade non-serrated - outer blade serrated

Ø 34 mm L 130 mm
Ø 4,2 mm L 130 mm
Ø 5,5 mm L 130 mm

LINVATEC

AU34AFR

AU42AFR

AU55AFR

ARTHREX

BU34AFR

BU42AFR

BU55AFR

DYONICS

CU34AFR

CU42AFR

CU55AFR

STRYKER

DU34AFR

DU42AFR

DU55AFR

AESCULAP

EU34AFR

EU42AFR

EU55AFR

FULL RADIUS RESECTOR
innen ungezahnt - außen ungezahnt
inner blade non-serrated - outer blade non-serrated

Ø 34 mm L 130 mm
Ø 4,2 mm L 130 mm
Ø 5,5 mm L 130 mm

LINVATEC

AU34FRR

AU42FRR

AU55FRR

ARTHREX

BU34FRR

BU42FRR

BU55FRR

DYONICS

CU34FRR

CU42FRR

CU55FRR

STRYKER

DU34FRR

DU42FRR

DU55FRR

AESCULAP

EU34FRR

EU42FRR

EU55FRR

INCISOR
innen gezahnt - außen gezahnt
inner blade serrated - outer blade serrated

Ø 34 mm L 130 mm
Ø 4,2 mm L 130 mm
Ø 5,5 mm L 130 mm

LINVATEC

AU34INC

AU42INC

AU55INC

ARTHREX

BU34INC

BU42INC

BU55INC

DYONICS

CU34INC

CU42INC

CU55INC

STRYKER

DU34INC

DU42INC

DU55INC

AESCULAP

EU34INC

EU42INC

EU55INC

ACROMIONIZER
ovaler Fräser
oval bur

Ø 4,2 mm L 130 mm
Ø 5,5 mm L 130 mm

LINVATEC

AU42ACR

AU55ACR

ARTHREX

BU42ACR

BU55ACR

DYONICS

CU42ACR

CU55ACR

STRYKER

DU42ACR

DU55ACR

AESCULAP

EU42ACR

EU55ACR

ROUND BUR
runder Fräser
round bur

Ø 4,2 mm L 130 mm
Ø 5,5 mm L 130 mm

LINVATEC

AU42ROB

AU55ROB

ARTHREX

BU42ROB

BU55ROB

DYONICS

CU42ROB

CU55ROB

STRYKER

DU42ROB

DU55ROB

AESCULAP

EU42ROB

EU55ROB

Ø = (Durchmesser Außenrohr; diameter outer tube)
L = (Nutzlänge; effective length)

click me



www.biq-shop.de

**K-wires
Configurator**



**Shaverblades
Configurator**



Alle hier aufgeführten Firmennamen und Handelsmarken sind immer Eigentum fremder Firmen – BiQ medical ist nicht affiliert mit all diesen Firmen! Alle Bilder und Abbildungen können vom Original abweichen. Änderungen bezüglich Inhalt und Darstellung vorbehalten! All company names and trademarks listed here are always the property of third-party companies – BiQ medical is not affiliated with any of these companies! All images and illustrations may differ from the original. Subject to change without notice in terms of content and presentation!

Shopping 24/7: www.biq-shop.de (including configurator for K-wires and Shaverblades).