



Micro saws

Mikrosägen mit Silikongriff

Griffsicherheit in trockener und feuchter Umgebung

Micro saws with silicone handle

Secure grip in dry and humid conditions

Micro sierras con mango de silicona

Seguridad de agarre en ambientes secos y húmedos

Micro scies avec poignée en silicone

Bonne prise en main dans un environnement sec ou humide

Micro seghe con impugnatura in silicone

Presu sicura in ambienti asciutti e umidi

Mikrosägen mit Silikongriff
Micro saws with silicone handle
Micro sierras con mango de silicona
Micro scies avec poignée en silicone
Micro seghe con impugnatura in silicone


Oszillierende Mikrosäge mit Schlüssel
 Micro oscillating saw with wrench
 Micro sierra oscilante con llave
 Micro scie oscillante avec clé
 Micro sega oscillante con chiave



INTRA		Sachse	68.24.01	12,5 cm - 5"		15.000 U/min

Sägeblätter
 Saw blades
 Hojas de sierra
 Lames de scie
 Lame da sega

Nutzlänge
 Working length
 Largo útil
 Longueur utile
 Lunghezza utile




	6 mm	10 mm	15 mm	6 mm
1:1				
	0,3 mm	0,3 mm	0,3 mm	0,15 mm
	a = 16 mm	a = 16 mm	a = 16 mm	a = 16 mm
P	68.24.24	68.24.25	68.24.26	68.24.30
	0,3 mm	0,3 mm	0,3 mm	0,15 mm
	a = 30 mm	a = 30 mm	a = 30 mm	a = 30 mm
P	68.24.27	68.24.28	68.24.29	68.24.31

Mikro-Stichsäge mit Schlüssel
 Micro compass saw with wrench
 Micro sierra de vaivén con llave
 Micro scie à guichet avec clé
 Micro sega con chiave



INTRA		Sachse	68.24.11	12,5 cm - 5"		15.000 U/min

Sägeblätter
 Saw blades
 Hojas de sierra
 Lames de scie
 Lame da sega



1:1					
	0,25 mm	0,25 mm	0,25 mm	0,25 mm	0,25 mm
	a = 11 mm	a = 14 mm	a = 18 mm	a = 22 mm	a = 26 mm
P	68.24.12	68.24.13	68.24.14	68.24.15	68.24.16



Maximale Motorendrehzahl = Oszillationen/Minute
 Maximum speed = oscillations/minute
 Revoluciones máximas del motor = oscilaciones/minuto
 Nombre de tours maximum = oscillations/minute
 Massimo numero di giri del motore = oscillazioni/minuto



Spülrohr, integriert
 Irrigation tube, integrated
 Tubo de irrigación, integrado
 Canule de lavage, intégrée
 Tubo di lavaggio, integrato



Spülrohr, abnehmbar
 Irrigation tube, detachable
 Tubo de irrigación, desmontable
 Canule de lavage, démontable
 Tubo di lavaggio, smontabile



Mikrosägen mit Silikongriff
Micro saws with silicone handle
Micro sierras con mango de silicona
Micro scies avec poignée en silicone
Micro seghe con impugnatura in silicone

Mikro-Sagittal-Säge mit Schlüssel
Micro sagittal saw with wrench
Micro sierra sagital con llave
Micro scie sagittale avec clé
Micro sega sagittale con chiave



Sachse

INTRA 


68.24.03

12,5 cm - 5"



20.000 U/min

Sägeblätter
Saw blades
Hojas de sierra
Lames de scie
Lame da sega



					
68.24.70	68.24.71	68.24.72	68.24.73	68.24.74	68.24.75

Osseoskalpell mit Schlüssel
Osseoscalpel with wrench
Osseoscalpel con llave
Osséoscalpel avec clé
Osseoscalpello con chiave



Sachse

INTRA 


68.24.06

12,5 cm - 5"




20.000 U/min

Sägeblätter
Saw blades
Hojas de sierra
Lames de scie
Lame da sega



1:1



		
68.24.35	68.24.36	68.24.37

Zubehör
Accessories
Accesorios
Accessoires
Accessori



68.24.19

Überwurfmutter für Spannfutter (allein)
Fixation nut for chuck (only)
Tuerca de racor para mandril de sujeción (sola)
Ecrou fileté pour mandrin de serrage (seule)
Dado per raccordi per mandrino (solo)



68.24.20

Gabelschlüssel (allein)
Wrench (only)
Llave (sola)
Clé (seule)
Chiave (sola)



68.24.22

Ersatz-Spühlrohr
Spare irrigation tube
Tubo para irrigación de recambio
Canule de lavage de rechange
Tubo di lavaggio di ricambio



Spühlrohr, abnehmbar
Irrigation tube, detachable
Tubo de irrigación, desmontable
Canule de lavage, démontable
Tubo di lavaggio, smontabile

Medicon eG

Gänsäcker 15 · D-78532 Tuttlingen · Germany · Tel.: +49 (0) 74 62 / 20 09-0
Fax: +49 (0) 74 62 / 20 09-50 · E-Mail: sales@medicon.de · Internet: www.medicon.de

medicon[®]

CMF

